


## Fast and Easy Tester Available for 80A09 TekProtect Protection Devices

### TF-TekProtect “Health Checker” Verifies EOS/ESD Protection is working

Verifying the TekProtect device in a customer’s 80A09 is still working properly is important to overall EOS/ESD protection for the unit. Manual procedures for verifying the TekProtect devices are included in Service Manuals for the products.



*TF-TekProtect ESD Protection Tester*

And a convenient “health checker” accessory is now available to quickly check the customer’s unit. The TF-TekProtect Health Checker is a small tester powered by USB connector, which connects to the 80A09. A single button-press tests the instrument’s TekProtect device in a few seconds and gives a Green (Pass) or Red (Fail) indication. Green means your instrument is protected; Red means your unit needs to have its TekProtect device replaced.